



Email : sales@laboratoryinstrumentindia.com

Product Name :
Portable Digital Press 56KN Capacity

Product Code :
LBNY-0005-1220007



Description :

The load is applied by a hand pump, and is measured by a high precision electric load cell with a digital display unit range 0-56 KN. Portable Digital Press 56KN Capacity Used for compression tests on small cylinder specimens and core samples up to dia. 60x100 mm. Providing:

0,001 KN resolution

65.000 divisions

Hysteresis: 0,03%

Repeatability: 0,02%

Linearity: 0,05%.

Technical Specification :

Complete with wooden carrying case, accessories.

The compression platens have dia. 65 mm., the upper one has a spherical seat and the vertical daylight is 110 mm.

Dimensions 370x320x710 mm.

Weight 25 kg.



Laboratory Instrument India

