



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Canin Instrument

**Product Code :**  
LBNY-0005-12200010



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

The large display just 9 keys for simple operation using menu technique and intelligent memory render Canin a unique instrument worldwide. 120.000 measurements can be stored in the intelligent memory and called up with the cursor keys. Canin corrosion analysing instrument, for the non-destructive detection of corrosion in the reinforcement bars of concrete building elements. Standard supply includes one bar electrode RS 232 integrated software for printer, cables, copper sulphate carrying case. A measuring surface of more than 4000 sq.mt.can be managed with the large memory.

**Technical Specification :**

Dimensions: 300x330x100 mm  
Weight: 5 Kg



**Laboratory Instrument India**