




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Ram Rapid Analysis Machine

**Product Code :**  
LBNY-0005-12200012



**Description :**

Fully automatic procedure with quick and accurate test results (max. errors within 5 kg./cubicmetre). Ram Rapid Analysis Machine Used for the determination of cement content in fresh concrete, coarse and fine aggregate, fly ash and GGBF slag content can also be determined. The complete test takes approx. 10 minutes. The connection to water net for approx. 80 liters each test is required.

**Technical Specification :**

Supplied complete.  
The weight of the test sample is 8 kg.  
Power supply: 230 V 1ph 50 Hz  
Dimensions: 780x660x1500 mm  
Weight: 150 kg



**Laboratory Instrument India**