




Email : sales@laboratoryinstrumentindia.com

Product Name :
Determination Of Toughness First Crack Strength And
Ductility Of Fibre Reinforced Concrete

Product Code :
LBNY-0005-12200017



Description :

The equipment can be used only with the flexure Servo-Plus Evolution testing machines models:
Flexure high stiffness frame 200kN capacity
Flexure open-side machine, 150kN capacity.

Technical Specification :

Determination Of Toughness First Crack Strength And Ductility Of Fibre Reinforced Concrete



Laboratory Instrument India