



Email : sales@laboratoryinstrumentindia.com

Product Name :
Slump Cone Galvanized Complete Set

Product Code :
LBNY-0005-12200018



Description :

Slump Cone complete set including:
Base plate, galvanized steel
Stainless steel rule, 300 mm long
Slump Cone, galvanized steel
Tamping rod, galvanized steel, dia. 16 x 600 mm
Fine wire brush
Aluminium scoop, 500 cc capacity

Technical Specification :

Weight: 5 kg.



Laboratory Instrument India