



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Slump Cone Plastic

**Product Code :**  
LBNY-0005-12200019



**Description :**

Slump Cone Plastic complete set including:  
Tamping rod, galvanized steel dia. 16 x 600 mm  
Base plate, galvanized steel  
Slump Cone, plastic Weight: 750 g  
Aluminium scoop, 500 cc capacity  
Stainless steel rule, 300 mm long  
Fine wire brush

**Technical Specification :**

Weight: 5 kg



**Laboratory Instrument India**