



Email : sales@laboratoryinstrumentindia.com

Product Name :
Slump Cone Plastic

Product Code :
LBNY-0005-12200019



Description :

Slump Cone Plastic complete set including:
Tamping rod, galvanised steel dia. 16 x 600 mm
Base plate, galvanised steel
Slump Cone, plastic Weight: 750 g
Aluminium scoop, 500 cc capacity
Stainless steel rule, 300 mm long
Fine wire brush

Technical Specification :

Weight: 5 kg



Laboratory Instrument India