



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
V Funnel

**Product Code :**  
LBNY-0005-12200020



**Description :**

The upper edge of the funnel is smooth and reinforced, and the outflow orifice is equipped of an open able seal valve. V Funnel stainless steel made, having 10 liters capacity stand mounted.

**Technical Specification :**

Dimensions: 510x400x1000 mm  
Weight: 40 kg



**Laboratory Instrument India**