



Email : sales@laboratoryinstrumentindia.com

Product Name :
Digital Concrete Test Hammer

Product Code :
LBNY-0005-12200021



Description :

The unit consists of the standard mechanical but equipped with an electronic transducer that measures the rebound values and supplies automatically the results on a graphic display. This Digital Concrete Test Hammer microprocessor operated, entirely designed and manufactured by with advanced technology, performs basic concrete testing with continuous automatic recording of all parameters in accordance with Specifications, register and process data and then transfer them to a PC.

Technical Specification :

Digital Concrete Test Hammer



Laboratory Instrument India