




Email : sales@laboratoryinstrumentindia.com

Product Name :
High Performance Ultrasonic Tester

Product Code :
LBNY-0005-12200027



Description :

High Performance Ultrasonic Tester is an instrument using the most modern technologies, it has a 14 VGA colour touch screen, 64 MB, Compact Flash interface, SD card, USB working system Windows CE with the possibility to manage, EXCEL, WORD, PPT files etc. Ultrasonic pulse velocity tester high performance with microprocessor for combined ultrasonic and rebound hammer data acquisition and processing.

Technical Specification :

High Performance Ultrasonic Tester



Laboratory Instrument India