



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Kelly Ball Apparatus

**Product Code :**  
LBNY-0005-12200030



**Description :**

The ball is lowered into the concrete and the penetration measured. Consisting of a hemispherical ally ended cylinder with guiding frame and a handle graduated in inch, it is used to determine the workability of fresh concrete.

Cadmium plated for rust protection.

It can be used on site or in laboratory.

**Technical Specification :**

Weight: 15 Kg



**Laboratory Instrument India**