



Email : sales@laboratoryinstrumentindia.com

Product Name :
Micro Coring Equipment

Product Code :
LBNY-0005-12200031



Description :

Micro Coring system is superiorly valid and reliable if combined with ultrasonic tester and concrete hammer. Micro-core extraction is easy, and requires the presence of one operator only. The extraction of a Micro Coring sample from a concrete structure or masonry is an extremely valid non-destructive method as it allows analysis and accurate evaluations of the manufacture without causing any damages to the structure, considering the dimension of the hole that can be eventually clogged with mortar.

Technical Specification :

Micro Coring Equipment comprises:

Flanged guide assembly

Drilling mask

Suitable electric drill. 230V 1F 50Hz

2 Self-blocking pincers to fit the flanged guide assembly to the surface

Impregnated diamond bit for cores with $\text{Ø} 28 \times 200 \text{ mm}$

Impregnated diamond bit for cores with $\text{Ø} 28 \times 100 \text{ mm}$

Set of accessories comprising: anchors, bits, wrenches, screws.

Carrying case.

Dimensions: 550x400x200 mm approx.

Weight: 10 Kg.



Laboratory Instrument India