




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
L Box Stainless

**Product Code :**  
LBNY-0005-12200032



**Description :**

L Box Stainless steel made consisting of:  
Container with inside rigid surfaces,  
Two marks far 200 and 400 mm are engraved on the horizontal bottom position.  
Obstacle for med by 3 vertical bars having 16 mm dia. with free light of 38 mm,  
gate in guillotine form.

**Technical Specification :**

Dimensions:700 x 300 x 650 mm  
Weight :40 kg



**Laboratory Instrument India**