



Email : sales@laboratoryinstrumentindia.com

Product Name :
Workshop Respirator

Product Code :
LBNY-0004-740022



Description :

Workshop Respirator

Technical Specification :

Respirator

Double filter tank and breathing valve design

Soft Material

Material - Made of soft silicon gel and active carbon, healthy, durable . It can protect you from damage of little dusts, mists, metallic fumes, spray paint

Healthy Material

Cartridge material: activated carbon

Great for filtering small dusts, mists, metallic fumes, spray paint, chemicals and small flying particles of low toxicity

Effectively filter poison gas and clean air

Made of healthy and clean material, without harmless

Easily Adjustable Strap

Easy to wear

Adjustable strap for head circumference of all sizes

Adjustable Belt: Comfortable and flexible headband and easy-to-adjust belt with neck clip for easy wear and tear and optimum pressure distribution.

High Quality Lens

The lenses are made of polyurethane, impact resistance.

Workshop Respirator, Workshop Respirator Equipments, Workshop Respirator Tools, Workshop Respirator Tool Kits, Workshop Respirator Manufacturers, Workshop Respirator Suppliers from India, China, Kenya



Laboratory Instrument India