



Email : sales@laboratoryinstrumentindia.com

Product Name :
Abrasion Tester Boehme

Product Code :
LBNY-0005-12200035



Description :

Perform the test for 16 cycles composed of 22 turn each, calculating at the end a worn as a average loss in volume and weight. Before doing a test, establish the specimens bulk density by measuring weight and thickness. The test is performed by positioning a specimen to be verified in a Abrasion Tester Boehme apparatus on the test track on which has been spread normalized abrasive the grinding wheel its made rotate and the specimen submitted to the abrasive load of 294 N for a certain number of cycles.

Technical Specification :

Abrasion Tester Boehme instrument measures a volume loss in a specimen under abrasion test and its used in tests such as:

Slabs made of natural rocks

Natural stone slabs

Paving stones

Concrete slabs

Power supply: 230V 50Hz 1PH 800 W

Dimension: 1500 x 1000 xh 850 mm

Weight: 250 kg

Laboratory Instrument India

