



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pull Out Test Apparatus

Product Code :
LBNY-0005-12200038



Description :

The standard equipment to Specification and comprises hydraulic extraction unit 100 kN capacity with pump, precision manometer 0-100kN, bearing ring, 10 steel discs 25 mm dia. carrying cases. Used to evaluate the concrete resistance as per the strength applied to extract a disc embedded into concrete.

Technical Specification :

Weight: 18 kg



Laboratory Instrument India