



Email : sales@laboratoryinstrumentindia.com

Product Name :
Vebe Consistometer

Product Code :
LBNY-0005-12200043



Description :

The Vebe Consistometer method is based on the same principle of the simple slump cone test method, but it has the advantage of a mechanized action, for the determination of the workability of concrete. After removing the slump cone, the concrete undergoes a vibration to determine its slump. Supplied complete.

Technical Specification :

Power supply: 230 V 1ph 50 Hz 250 W
Dimensions: 260x380x700 mm
Weight: 90 kg.



Laboratory Instrument India