



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Analogic Strain Gauge 300 MM Length

**Product Code :**  
LBNY-0005-12200044



**Description :**

Analogic Strain Gauge 300 MM Length instrument can also be used for other structures like steel and wood. Strain gauge, 300 mm. measuring length. Different models are available with analogic or digital gauge, 100, 200, 300 mm measuring length, depending on the standard length to be measured. Analogic Strain Gauge 300 MM Length Used to determine the strain (length changes) in concrete specimens and structures, rock strata, different parts of a structure, in remote areas and under adverse conditions, using a single instrument.

**Technical Specification :**

Analogic Strain Gauge 300 MM Length



**Laboratory Instrument India**