



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Drum Type Mixer

**Product Code :**  
LBNY-0005-12200045



**Description :**

Suitable for field mixes of low medium strength concrete.

Drum volume: 125 liters

Power supply: 230V 1ph 50 Hz - 0,3 HP

Yield: 75 liters of concrete

**Technical Specification :**

Dimensions: 720x1320x1280 mm

Weight: 60 Kg.



**Laboratory Instrument India**