



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Eight Channels Control Unit Touch Screen Semi Automatic

**Product Code :**  
LBNY-0005-12200048



**Description :**

The new PC based and touch screen control unit developed by department to offer a new generation of control units for the management of semi-automatic and fully automatic compression, flexural, and splitting tensile tests on concrete and cement specimens.

**Main Features:**

- Modular structure
- Pc based control unit
- No limits memory
- Touch screen display with icons Internet and Intranet connections for remote technical assistance.

**Technical Specification :**

Eight Channels Control Unit Touch Screen Semi Automatic



**Laboratory Instrument India**