



Email : sales@laboratoryinstrumentindia.com

Product Name :
Eight Channels Control Unit Automatic

Product Code :
LBNY-0005-12200053



Description :

Designed with the latest technology, an innovative PC-like Touch Screen system, employed to control and manage all sorts of automatic testing machines. Multi piston electric pump with variable flow driven by a microprocessor. To update or complete your compression and flexural testing machine on concrete and mortar. Eight channels controlled unit for a fully automatic execution of the test.
Data acquisition system.
The machine comprises.

Technical Specification :

Eight Channels Control Unit Automatic



Laboratory Instrument India