



Email : sales@laboratoryinstrumentindia.com

Product Name :
Workshop Welding Apron

Product Code :
LBNY-0004-740024



Description :

Workshop Welding Apron

Technical Specification :

Welding Apron
Standard for welding
Classification: Sleeveless Apron
Material: Leather
Size: Free size
Pattern: Plain
Disposable: Disposable
Gender: Unisex.

Workshop Welding Apron, Workshop Welding Apron Equipments, Workshop Welding Apron Tools, Workshop Welding Apron Tool Kits, Workshop Welding Apron Manufacturers, Workshop Welding Apron Suppliers from India, China, Kenya



Laboratory Instrument India

