



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Pan Type Mixer Yield 55 Liters

**Product Code :**  
LBNY-0005-12200055



**Description :**

Yield per mixture lt. 55  
Motor power 1,1 KW  
Cylinder dimensions (dia x h) 70x30 cm  
Cylinder capacity lt. 100

**Technical Specification :**

Dimensions (dia x h) 71x115 cm  
Weight kg 115.



**Laboratory Instrument India**