




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Flexural Testing Machine 150kN Digital Model Open Side  
Frame

**Product Code :**  
LBNY-0005-12200058



**Description :**

Flexural open sided frame 150 KN capacity complete with pressure transducer used in conjunction with a compression machine equipped with digital control units.

**Technical Specification :**

Flexural Testing Machine 150kN Digital Model Open Side Frame



**Laboratory Instrument India**