




Email : sales@laboratoryinstrumentindia.com

Product Name :
Curing Tank 1000 Liters Capacity

Product Code :
LBNY-0005-12200061



Description :

Made from steel sheet, zinc coated to prevent it from corrosion. The tank can accommodate up to 64 cubes 150mm side, or up to 48 cubes 200mm side. Complete with base rack and stopper for an easy water discharge. Supplied without thermostat heating system, to be ordered separately.

Technical Specification :

Inside dimensions: 1500x750x750 mm
Weight: 120 kg.



Laboratory Instrument India