



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Split Cylinder Mould Dia. 160x320 MM

**Product Code :**  
LBNY-0005-12200062



**Description :**

Foreseen of lateral hinges for total opening and fast clamping system with inbuilt revolving screw. Steel made, galvanized finishing against corrosion. Complete with base they are easy to use with practical and fast demoulding recommended for field use.

**Technical Specification :**

Dimensions: dia. 160x320 mm  
Weight: 11 kg



**Laboratory Instrument India**