



Email : sales@laboratoryinstrumentindia.com

Product Name :
Electro Hydraulic Loading And Control System

Product Code :
LBNY-0005-12200063



Description :

Electro Hydraulic loading and control system consisting of:
Upper attachment for steel frame cross-mean coupling.
Spherical seat fixed to the ram for an uniform loading.
Double action alloy steel ram + cylinder. The ram is ground.

Technical Specification :

Max. ram travel: 400 mm
Max. load capacity: 1000 kN.



Laboratory Instrument India