



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Electro Hydraulic Loading And Control System

**Product Code :**  
LBNY-0005-12200063



**Description :**

Electro Hydraulic loading and control system consisting of:  
Upper attachment for steel frame cross-mean coupling.  
Spherical seat fixed to the ram for an uniform loading.  
Double action alloy steel ram + cylinder. The ram is ground.

**Technical Specification :**

Max. ram travel: 400 mm  
Max. load capacity: 1000 kN.



**Laboratory Instrument India**