



Email : sales@laboratoryinstrumentindia.com

Product Name :
Workshop Gauntlet Cuff Gloves

Product Code :
LBNY-0004-740025



Description :

Workshop Gauntlet Cuff Gloves

Technical Specification :

Gauntlet Cuff Gloves
Standard for Welding
Size: 14 Inch
Feature: Breathable, Heat Resistant, Durable
Style: Leather Gloves

Workshop Gauntlet Cuff Gloves, Workshop Gauntlet Cuff Gloves Equipments, Workshop Gauntlet Cuff Gloves Tools, Workshop Gauntlet Cuff Gloves Tool Kits, Workshop Gauntlet Cuff Gloves Manufacturers, Workshop Gauntlet Cuff Gloves Suppliers from India, China, Kenya



Laboratory Instrument India