



Email : sales@laboratoryinstrumentindia.com

Product Name :
Accelerated Concrete Curing Tank

Product Code :
LBNY-0005-12200065



Description :

Accelerated Concrete Curing Tank can hold up to 16 cubic 150 mm.side specimens, or 16 cylindrical dia.150 mm.specimens, or 8 cubic 200 mm.side specimens.

Accelerated Concrete Curing Tank comprises a fully insulated double wall tank with cover, inside all from stainless steel, outside from steel painted sheet with an intermediate layer of insulating mineral wool. Accelerated Concrete Curing Tank has been designed for accelerated concrete strength curing.

The test consists essentially in curing the concrete specimens with water heated by 3 electric elements of 1500W each.

The separate control panel is provided with a thermo regulator, timer, pilot lights, main switch. Temperature range: from ambient to 100 °C.

Technical Specification :

Power supply: 220-240V 1ph 50/60 Hz 4500W

Inside dimensions: 910x660x680 mm

Overall dimensions: 970x720x900 mm

Weight: 130 kg.

Laboratory Instrument India

