



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pan Type Mixer 56 Liters

Product Code :
LBNY-0005-12200068



Description :

The pan is easily removable by means of a trolley. The blades are hardened against wear. This multiflow mixer absorbs fewer air during mixing, requires shorter mixing time and grants a perfect homogeneity in mixtures having a low water cement ratio.

Technical Specification :

Power supply: 220-240 V 1ph 50 Hz 2 Hp
Mixing pan: 640 mm.dia.x 330 mm deep
Weight: 250 Kg.



Laboratory Instrument India