



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Pan Type Mixer 56 Liters

**Product Code :**  
LBNY-0005-12200068



**Description :**

The pan is easily removable by means of a trolley. The blades are hardened against wear. This multiflow mixer absorbs fewer air during mixing, requires shorter mixing time and grants a perfect homogeneity in mixtures having a low water cement ratio.

**Technical Specification :**

Power supply: 220-240 V 1ph 50 Hz 2 Hp  
Mixing pan: 640 mm.dia.x 330 mm deep  
Weight: 250 Kg.



**Laboratory Instrument India**