



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Two Channels Control Unit Semi Automatic

**Product Code :**  
LBNY-0005-12200069



**Description :**

Two channels computerised graphic display system to control and manage all sorts of semi-automatic testing machines, processing, printing and saving the test data and certificates, for acquisition, display, with software for remote control from PC.

**Technical Specification :**

Automatic elaboration of the specific resistance value.  
Permanent file up to 1000 tests and file of 100 different types of specimens.  
Simple and immediate set up of the parameters and test execution, menu driven. The use does not require specialised staff.  
Two analogue-digital channels connectable to two different compression/flexure frames.  
Class: 0,5% starting from 10% of maximum value, on request from 1% of maximum value.  
Continue load display.  
Breaking load detection.  
Selectable measuring force: KN, lb  
Graphic display with high resolution: 192x64 pixels.



**Laboratory Instrument India**

