



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Curing Tank 650 Liters Capacity Heavy Plastic

**Product Code :**  
LBNY-0005-12200070



**Description :**

Made from extremely robust and stable polyethylene, complete with base rack. Curing tank 650 liters capacity, heavy plastic. Supplied without thermostat heating system to be ordered separately.

**Technical Specification :**

Inside dimensions: 1040 x 1040 x 605 mm  
Weight: 60 kg.



**Laboratory Instrument India**