



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Workshop TIG Welding Gloves

**Product Code :**  
LBNY-0004-740026



**Description :**

Workshop TIG Welding Gloves

**Technical Specification :**

TIG Welding Gloves  
Standard for welding  
Material: Leather  
Feature: Abrasion Resistant, Heat Resistant  
Size: 16 Inch (40CM)

Workshop TIG Welding Gloves, Workshop TIG Welding Gloves Equipments, Workshop TIG Welding Gloves Tools, Workshop TIG Welding Gloves Tool Kits, Workshop TIG Welding Gloves Manufacturers, Workshop TIG Welding Gloves Suppliers from India, China, Kenya



**Laboratory Instrument India**