



Email : sales@laboratoryinstrumentindia.com

Product Name :
Workshop TIG Welding Gloves

Product Code :
LBNY-0004-740026



Description :

Workshop TIG Welding Gloves

Technical Specification :

TIG Welding Gloves
Standard for welding
Material: Leather
Feature: Abrasion Resistant, Heat Resistant
Size: 16 Inch (40CM)

Workshop TIG Welding Gloves, Workshop TIG Welding Gloves Equipments, Workshop TIG Welding Gloves Tools, Workshop TIG Welding Gloves Tool Kits, Workshop TIG Welding Gloves Manufacturers, Workshop TIG Welding Gloves Suppliers from India, China, Kenya



Laboratory Instrument India