



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Specimen Bench Cutting Machine

**Product Code :**  
LBNY-0005-12200074



**Description :**

The machine accepts blades up to dia. 350 mm  
Blade rotation speed: 3900 rpm  
Shear capacity: 120 mm  
Power supply: 230 V 1ph 50 Hz 2000 W  
Supplied complete with abrasive blade dia. 350 mm

**Technical Specification :**

Dimensions: 560x460x390 mm  
Weight: 20 Kg.



**Laboratory Instrument India**