



Email : sales@laboratoryinstrumentindia.com

Product Name :
Specimen Bench Cutting Machine

Product Code :
LBNY-0005-12200074



Description :

The machine accepts blades up to dia. 350 mm
Blade rotation speed: 3900 rpm
Shear capacity: 120 mm
Power supply: 230 V 1ph 50 Hz 2000 W
Supplied complete with abrasive blade dia. 350 mm

Technical Specification :

Dimensions: 560x460x390 mm
Weight: 20 Kg.



Laboratory Instrument India