



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Concrete Flexural Machine 150 KN Motorized

**Product Code :**  
LBNY-0005-12200075



**Description :**

Concrete flexural machine 150 KN motorized with control unit to perform flexural tests on concrete beam specimens max. dimensions 150x150x750 (600) mm.

**Technical Specification :**

Concrete Flexural Machine 150 KN Motorized



**Laboratory Instrument India**