



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Concrete Flexural Machine 150 KN Open Sided Frame  
Motorized Fully Automatic Model

**Product Code :**  
LBNY-0005-12200076



**Description :**

Concrete flexural machine 150 KN open sided frame motorized with control unit fully automatic model to perform flexural tests on concrete beam specimens max. dimensions 200x200x800mm, kerbs, tiles, slabs, masonry units, flat blocks, flagstones, and any type of material having max. size 600x250 mm (lower rollers max. length 1325 mm).

**Technical Specification :**

Concrete Flexural Machine 150 KN Open Sided Frame Motorized Fully Automatic Model



**Laboratory Instrument India**