




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Concrete Flexural Machine 200 KN High Stiffness Semi Automatic Model

**Product Code :**  
LBNY-0005-12200078



**Description :**

Concrete flexural machine 200 kN, high stiffness with control unit (semi-automatic model), to perform flexural tests on concrete beam specimens max. dimensions 200x200x800mm, kerbs, tiles, slabs, flat blocks, flagstones, masonry units, and any type of material having max. size 600x250 mm (lower rollers max. length 1325 mm).

**Technical Specification :**

Concrete Flexural Machine 200 KN High Stiffness Semi Automatic Model



**Laboratory Instrument India**