



Email : sales@laboratoryinstrumentindia.com

Product Name :
Concrete Flexural Machine 150 KN Motorized Fully Automatic Model

Product Code :
LBNY-0005-12200080



Description :

Concrete flexural machine 150 KN motorized with Servo-Plus evolution unit fully automatic model to perform flexural tests on concrete beam specimens max. dimensions 150x150x750 (600) mm.

Technical Specification :

Concrete Flexural Machine 150 KN Motorized Fully Automatic Model



Laboratory Instrument India