



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Concrete Flexural Machine 150 KN Motorized Fully Automatic Model

**Product Code :**  
LBNY-0005-12200080



**Description :**

Concrete flexural machine 150 KN motorized with Servo-Plus evolution unit fully automatic model to perform flexural tests on concrete beam specimens max. dimensions 150x150x750 (600) mm.

**Technical Specification :**

Concrete Flexural Machine 150 KN Motorized Fully Automatic Model



**Laboratory Instrument India**