



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Specimen Cutting Machine Sliding Supports Model

**Product Code :**  
LBNY-0005-1220081



**Description :**

Specimen cutting machine with sliding supports  
Useful cutting height: 130 mm  
The machine accepts blades up to dia. 400 mm  
Blade rotation speed: 2800 rpm  
Supplied without blade

**Technical Specification :**

Power supply: 230 V 1ph 50 Hz 3 Hp  
Dimensions of the sliding table: 460X400 mm  
Dimensions: 1185x660x1400 mm  
Weight: 79 kg.



**Laboratory Instrument India**