



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Concrete Compression Machine 5000 KN Motorized

**Product Code :**  
LBNY-0005-12200082



**Description :**

Compression testing machine 5000 KN motorized with control unit (fully automatic model), to test cubes up to 300 mm side and cylinders up to dia. 250x500 mm.

**Technical Specification :**

Concrete Compression Machine 5000 KN Motorized



**Laboratory Instrument India**