



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Specimen Cutting Machine Blade Max. Dia. 600 MM

**Product Code :**  
LBNY-0005-12200083



**Description :**

Specimen Cutting Machine can accept blade having max. dia. 600 mm.  
Power supply: 400 V 3ph 50 Hz 5,5 Hp

**Technical Specification :**

Specimen Cutting Machine Blade Max. Dia. 600 MM



**Laboratory Instrument India**