



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Compression Flexural Frame On Mortar Specimens

**Product Code :**  
LBNY-0005-12200084



**Description :**

Compression Flexural frame on mortar specimens, dual range:  
0-15 kN for flexure tests complete with two pressure transducers used along with a digital concrete compression machine  
0-250 kN for compression tests.

**Technical Specification :**

Compression Flexural Frame On Mortar Specimens



**Laboratory Instrument India**