



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Concrete Compression Machine 5000 KN High Stability

**Product Code :**  
LBNY-0005-12200086



**Description :**

Concrete compression machine 5000 KN high stability with control unit the ideal solution for central and research laboratories to test high strength specimens explosive samples rock and ceramic samples.

**Technical Specification :**

Concrete Compression Machine 5000 KN High Stability



**Laboratory Instrument India**