




Email : sales@laboratoryinstrumentindia.com

Product Name :
Concrete Flexural Machine 200 KN High Stiffness Fully Automatic Model

Product Code :
LBNY-0005-12200087



Description :

Concrete flexural machine 200 KN high stiffness, with control unit (fully automatic model) to perform flexural tests on concrete beam specimens max. dimensions 200x200x800mm, kerbs, tiles, slabs, masonry units, flat blocks, flagstones, and any type of material having max. size 600x250 mm (lower rollers max. length 1325 mm).

Technical Specification :

Concrete Flexural Machine 200 KN High Stiffness Fully Automatic Model



Laboratory Instrument India