




Email : sales@laboratoryinstrumentindia.com

Product Name :
Concrete Flexural Machine 150 KN Open Sided Frame
Motorized Semi Automatic Model

Product Code :
LBNY-0005-12200088



Description :

Concrete flexural machine 150 kN open sided frame, motorized with control unit (semi-automatic model) to perform flexural tests on concrete beam specimens max. dimensions 200x200x800mm, tiles, slabs, masonry units, flat blocks, flagstones, kerbs, and any type of material having max. size 600x250 mm (lower rollers max. length 1325 mm).

Technical Specification :

Concrete Flexural Machine 150 KN Open Sided Frame Motorized Semi Automatic Model



Laboratory Instrument India